




Raddec - New innovations

Phil Warwick



Pyrolyser Gen III




New front panel layout


All controls are now at the front

All electronics are housed within a sliding tray accessed at the front of the furnace

Pyrolyser Gen III draw



Straightforward maintenance and repair from the front of the furnace



Pyrolyser Gen IV


Improved system building and servicing by using an instrument drawer system

Enhanced electrical build to meet CSA standard

24 Volt DC components for enhanced safety

Single centrifugal fan with manifold to cool two zones

AIR, OXYGEN supplies - NITROGEN OPTION







Pyrolyser-6 Trio Generation^{IV}

- Incorporates a EUROTHERM HMI-PLC system
- E+PLC400 racks with modules
- 7" HMI (touch sensitive programming)
- Three EPC2000 O/T controllers
- Sample and Mid-zone rapid cool down facility for rapid cycling (fan+chimney)

The **Pyrolyser- Trio GenIV HMI-PLC** system :-

- 3 completely independent furnace zones
- Stores 20 editable heating programs;
- each program offers up to 15 segments.
- Ability to log and store combustion parameter data.
- The user can modify existing programs whilst the Pyrolyser is running a current program.
- A multi-level LOGIN system is available to manage users at different levels (Operator, Supervisor, Engineer)

The **Pyrolyser- Trio GenIV HMI-PLC** system stores 20 editable programs; each program offers up to 15 segments. The user can modify existing programs whilst the Pyrolyser is running a current program. A multi-level LOGIN system is available to manage users at different levels (Operator, Supervisor, Engineer)




The 'Data Entry' screen displays a table for entering data. At the top, it says 'Measuring Program:' and 'Analyst:'. The table has columns for 'Sample Name', 'Sample Type', 'Wet or Dry', and 'Sample Mass Grams' across six 'Worktube' columns (1 to 6). The 'Sample Name' and 'Sample Mass Grams' columns contain 'Test3' and '1.00' respectively. The 'Sample Type' and 'Wet or Dry' columns contain 'New Samp' and 'Wet'. There are also 'Copy' and 'Worktube 1' buttons at the bottom right.

	Worktube 1	Worktube 2	Worktube 3	Worktube 4	Worktube 5	Worktube 6
Sample Name	Test3	Test3	Test3	Test3	Test3	Test3
Sample Type	New Samp	New Samp	New Samp	New Samp	New Samp	New Samp
Catalyst Type	0.3% Pt-Al	0.3% Pt-Al	0.3% Pt-Al	0.3% Pt-Al	0.3% Pt-Al	0.3% Pt-Al
Wet or Dry	Wet	Wet	Wet	Wet	Wet	Wet
Sample Mass Grams	1.00	1.00	1.00	1.00	1.00	1.00

The 'Alarms' screen displays a table of active alarms. The table has columns for 'Active Date', 'Active Time', 'AlarmMessage', and 'AlarmStatus'. There are three active alarms listed.

	Active Date	Active Time	AlarmMessage	AlarmStatus
Active	06/07/2021	11:47:11	Sample Zone Over Temp	Active
Historical	06/07/2021	11:47:11	Mid Zone Over Temp	Active
Historical	06/07/2021	11:47:11	Catalyst Zone Over Temp	Active

Pyrolyser mini



Modular construction for installation in glovebox environments.


Readily accessible components for ease of maintenance.

Elements installed in a replaceable module.

Remotely located control unit

Multi-function capability.

Pyrolyser Mini trials





Samples loaded at 600°C then heated to 900°C.

Catalyst zone held at 800°C

Warwick P.E., Croudace I.W. & Burrell F.M. A compact, dual-zone vertical tube furnace for the determination of tritium and carbon-14 in decommissioning wastes. (2021) *Appl. Radiat. Isot.* **179**, 109995.

This work was supported by the Nuclear Decommissioning Authority under their Direct Research Portfolio [NS4510-500-003].

^3H evolution profiles




Extraction of ^3H and ^{14}C using a Pyrolyser Mini

Sample	H-3 activity (Bq/g)			C-14 activity (Bq/g)	
	Working value	Measured in 0.1M HNO_3	Measured in Carbontrap	Working value	Measured
Structural concrete (1.5 g)	4.8 ± 0.5	6.5 ± 0.8	6.9 ± 0.8	0.55 ± 0.07	0.59 ± 0.08
Irradiated concrete 1 (1.5 g)	37 ± 8	38 ± 4	46 ± 5	2.5 ± 0.3	2.8 ± 0.3
Irradiated concrete 2 (1.5 g)	1700 ± 200	1900 ± 200	2200 ± 300	7.7 ± 0.9	8 ± 1
Swarf metal 1 (1.5 g)	0.9 ± 0.6	7.6 ± 0.9	2.2 ± 0.3	7.0 ± 0.3	8 ± 1
Swarf metal 2 (1.5 g)	4 ± 7	3.4 ± 0.4	19 ± 2	7.9 ± 0.1	8 ± 1
OBT sediment sample (0.5 g)	11 ± 3	10 ± 1	10 ± 1	0.44 ± 0.04	0.6 ± 0.1

Note that the swarf metal contamination is highly heterogeneous


Accessories



Integrated vial bubbler

Pt-alumina catalyst

Introduction of 0.3% Pt-alumina catalyst



More robust supply of material.
0.5% Pt-alumina is being phased out.

Tested with ^3H -thymidine spiked milk.

Catalyst operating at 400°C.

Quantitative recovery of ^3H


No colouration in the bubbler.

Quench values comparable to those obtained using 0.5% Pt-alumina

Other improvements

- New standard worktube with a B34 entry cone to enable larger diameter boats (and hence sample) to enter.
- Provision of larger diameter worktubes (up to 42 mm OD to enable even larger sample to be loaded (e.g. 10 grams Fish, biota, foodstuffs)
- Development of a method to cleanly oxidise 10 grams fish in 7.5 hours (so 60 grams of fish in 1 working day)

Training



Online

- Main web site at
www.raddec.com
- Also, see technical videos on YouTube
(follow link from our website)

Acknowledgements

Pyrolyser Mini ^3H & ^{14}C extraction studies

Jonathan Cox (formerly NSG Environmental) and David Wickenden (MagnaX Ltd), for their oversight of this research

Ben Russell (National Physical Laboratory) for a contribution relating to the mobile laboratory deployment and

Yvonne Morris and the late David Loudon (Nuclear Decommissioning Authority) for their support of the study.

RADDEC

INTERNATIONAL



www.raddec.com