

# Getting the most out of your PT Results

- Why participate in proficiency tests?
- How to evaluate the results?
- How to implement the evaluation in your QA management system?




# Clever use of PT results can save time


---

- Quality Assurance Management Systems (such as ISO-17025, ISO-9001, ISO-15189 and others) impose strict requirements to:
  - Method validation
  - Method verification
  - Uncertainty budget
- General guidelines on method validation (such as Eurochem Guide 'Fitness for Purpose') specify a number of parameters:
  - Selectivity
  - Limit of Detection or Limit of Quantification
  - Linearity and Working Range
  - Trueness
  - Precision (Repeatability and Reproducibility)
  - Uncertainty
  - Ruggedness
- Clever use of PT results can provide information for a number of these parameters


# OBT PT #6: Large variation in the reported uncertainty




# Investigating the possible bias of your result



# Using your lab performance score (zèta, U, E<sub>n</sub>, ...)



# Effect of uncertainty on performance indicator $E_n$



# Getting more out of Proficiency Test results

- How to verify if your lab results ...
  - ... are correct
  - ... have a realistic estimation of uncertainty
  - ... are reproducible
  - ... are not biased




- An alternative method was proposed by Meijer et al. (\*) in 2002, based on the long term analytical performance of a laboratory in various proficiency tests (LTUM: Long Term Uncertainty Method). This method was applied by Matar er al (\*\*) in 2015.

(\*) Meijer P., de Maat M.P., Klufft C., Haverkate F. and van Houwelingen H.C.; "Long-term analytical performance of hemostasis field methods as assessed by evaluation of the results of an external quality assessment program for antithrombin."; Clin. Chem.; 2002; **48**: 1011-5.


(\*\*) Matar G., Poggi B., Mely R., Bon C., Chardon L., Chikh K., Renard A.C., Sota C., Eynard J.C., Cartier R. and Cohen R.; "Uncertainty in measurement for 43 biochemistry, immunoassay, and hemostasis rou)ne analytes ecaluated by a method using only external quality assessment data."; Clin. Chem. Lab. Med.; 2015; **53** (11); 1725-36.




# The 'perfect' lab: equality




# Effect of a constant bias



# Effect of a proportional bias



# The 'real' lab: combined (total) uncertainty



# The mathematics behind the analysis

$$TE = \sqrt{LTCV^2 + BT^2}$$

*'sum of squares'*

$$BT = \sqrt{BC^2 + BP^2}$$

$$BC = \sqrt{(\bar{y} - \bar{x})^2} = a$$


*~ intercept*

$$BP = \sqrt{\frac{(n-1)}{n} \times (b-1)^2 \times S_x^2}$$

*~ slope*

$$RE = \sqrt{\frac{(n-2)}{n} \times S_{\frac{y}{x}}^2}$$

$$S_{\frac{y}{x}} = \sqrt{\frac{1}{(n-2)} \times \left[ \sum (y - \bar{y})^2 - \frac{[\sum (x - \bar{x})(y - \bar{y})]^2}{\sum (x - \bar{x})^2} \right]}$$




Long term analytical error

Residual standard deviation (variability of the regression line)


$$LTCV = \frac{S_{\frac{y}{x}}}{\bar{x}} \times 100$$

Long term coefficient of variation (CV)

# Liquid Scintillation Counting: $^3\text{H}$ in OBT



# Liquid Scintillation Counting: $^{14}\text{C}$ in various matrices



## Conclusion: PT results provide excellent validation and verification data

---

- Easy to use method based on simple regression analysis
- Directly compare the performance of your lab with other labs
- Good estimate of long term stability (  $\Rightarrow$  Reproducibility)
- Good estimate of linearity (and possibly range)
- Good estimate of bias (both constant and proportional)
- Confirmation of calculated uncertainty budget
- Confirmation about trueness and precision
- Possible information about method ruggedness



# Questions?

---



## Copyright © SCK CEN

### PLEASE NOTE!

This presentation contains data, information and formats for dedicated use only and may not be communicated, copied, reproduced, distributed or cited without the explicit written permission of SCK CEN.

If this explicit written permission has been obtained, please reference the author, followed by 'by courtesy of SCK CEN'.

Any infringement to this rule is illegal and entitles to claim damages from the infringer, without prejudice to any other right in case of granting a patent or registration in the field of intellectual property.

### **SCK CEN**

Belgian Nuclear Research Centre

Foundation of Public Utility

Registered Office: Avenue Herrmann-Debrouxlaan 40 – BE-1160 BRUSSELS

Operational Office: Boeretang 200 – BE-2400 MOL